

The study of circadian rhythms and suprachiasmatic nucleus (SCN) neurons

How does studying jet lag in mice improve our health? A recent study led by Huiyan Li, from the National Center of Biomedical Analysis in Beijing, China is looking to understand how Chronobiology affects our wellbeing.

Circadian rhythms are physical, mental, and behavioral changes that follow a 24-hour cycle. This rhythm is governed by the master circadian pacemaker, suprachiasmatic nucleus (SCN) which can orchestrate the peripheral clocks in multiple tissues throughout the body. They have established primary cilia-mediated Sonic Hedgehog (SHH) signaling in the SCN as a novel interneuronal coupling mechanism and may lead to novel therapy of circadian disruption-linked diseases like high blood pressure, obesity and other metabolic disorders.

In vertebrate animals, the master clock is a group of about 20,000 nerve cells (neurons) that form a structure called SCN, which is in a part of the brain called the hypothalamus and receives direct input from the eyes. It is the central pacemaker of the circadian timing system and regulates most circadian rhythms in the body. Your body tries to align your sleep-wake cycle to cues from the environment, such as when it gets light or dark outside, when you eat, and when you are physically active.

Environmental circadian disruptions, such as acute jet-lag and long-term shift work, cause temporal unsynchronization between the internal circadian clock and external time cues, leading to physiological stress. Circadian disruption has been implicated in tumorigenesis and various psychiatric, neurological and metabolic diseases, including depression and diabetes. The study of circadian rhythms and suprachiasmatic nucleus (SCN) neurons is an active area of research, and new discoveries continue to expand our understanding of these complex processes.




Fig 1 - Primary Cilia in the SCN exhibit circadian dynamics
*Fig. 1. (A) Schematic diagram of the suprachiasmatic nucleus (SCN). (B) Representative three-dimensional reconstructed projection images of the SCN at 20 \times magnification of the two-photon imaging are shown. SCN slices were stained with anti-ACIII antibody (primary cilia marker, green). Scale bars, 50 μ m. (C) Representative images of primary cilia and expression of clock gene *Cry1* in the SCN during light/dark (LD) cycle. Scale bars, 100 μ m (main image) and 20 μ m (magnified region).*

Primary cilia, specialized hair-like structures are found on the surface of many cells, play crucial role in mammalian embryonic development. Tu et al. revealed that primary cilia in a subset of SCN neurons cilia exhibit pronounced circadian rhythmicity in abundance and length, and further identify primary cilia as a critical device for intercellular coupling to maintain the circadian clock in the SCN.




Fig 2 - Representative time-lapse images of the primary cilium in SCN neurons for 48 hours

Isolated live SCN neurons from postnatal mice were transduced with modified *baculovirus* encoding *mCherry* tagged, constitutively ciliary-localized protein, 5-hydroxytryptamine receptor 6 (5-HT6). The numbers on the images indicate the time in hours. Arrows indicate primary cilia. Scale bar, 5 μm . (J) Quantitative analysis of the cilium length in (I). Each dot represents one cell.

In vertebrates, primary cilia are required for transactivation of the Sonic Hedgehog (SHH) signalling. SHH pathways play a significant role in embryonic development and tissue patterning, primarily during early development. The research highlighted that cilia-mediated SHH signalling in the SCN is essential to control the central clock function in mice. They also found that a clinical drug, Vismodegib, blocking SHH signaling, could be used to regulate circadian rhythm, and render mice susceptible to rapid resetting by light. Thus, targeting cilia-mediated SHH signaling might be a potential therapeutic strategy for the treatment of human diseases related to circadian disruptions.

To test whether primary cilia-mediated SHH signaling is required for intercellular coupling among SCN neurons, the researchers used real-time luciferase luminescence imaging of SCN slices isolated from *Per2::Luciferase* (*Per2::Luc*) transgenic reporter mice to track *Per2* rhythmic expression in single cells *ex vivo*. In the experimental set-up bioluminescence imaging was used with identical culture conditions for SCN slices with a higher concentration of luciferin (1mM), which was added for imaging. The culture dish was sealed and placed on the stage of an inverted microscope with 10x objective lens in a dark room. See Fig 3. Images were acquired with three deep-cooled CCD cameras: Andor iKon-M 934, Raptor Eagle 47-10 and Raptor Eagle CCD 42-40. Images of 60 min exposure duration were collected continuously. The Eagle CCD42-40 camera came out as the winner with the optimum performance.

The team are well on their way to understanding the complexities of circadian rhythms and how they affect our wellbeing in different environments. They will look to develop therapeutic treatments.



Fig 3 - Experimental camera set-up

Raptor has been developing deep cooled CCD cameras for years. The Eagle family of cooled cameras enables longer integration times. Both 4MP and 1MP models are available offering -70°C and -90°C absolute cooling.



Fig 4 - Eagle camera

- **7-year vacuum guarantee** – Protection and integrity of the sensor
- **Extremely low dark current** – Deep cooled to greater than -110°C delta enables long exposure times
- **Back illuminated 4MP sensors from e2v** – Enables large field of view imaging
- **13.5 μm x 13.5 μm pixels (4MP)** – Enables ultra sharp image resolution
- **High QE: >90% @ 525nm and 50% @ 380nm & 720nm** – Optimum photon collection

Ref 1 - Rhythmic cilia changes support SCN neuron coherence in circadian clock- Hai-Qing Tu et al - National Center of Biomedical Analysis in Beijing, China.

doi: [https:// DOI: 10.1126/science.abm1962](https://doi.org/10.1126/science.abm1962)

Raptor UK (Headquarters)
T: +44(0)2828 270 141
E: sales@raptorphotronics.com
www.raptorphotronics.com

Raptor USA
T: +1 (877) 240-4836
E: sales@raptorphotronics.com
www.raptorphotronics.com

